

## DIAGNOSTIC WAX UP CHECK LIST

- High precision impression preferably in Polyvinyl or Polyether material
- Face Bow, taken parallel to the floor, patient in upright position. The use of a leveler is suggested. Mounting of the upper arch in the office is suggested to avoid any movements during shipping.
- Aesthetic photos (as per Eurodent protocol).
- Opposing impression or bubble free model poured in die stone.
- Bite relation(s) ensuring that CR and CO coincide. If edentulous areas are present, please adopt a base plate and a bite rim to ensure stability.
- If a therapeutic splint was used, please enclose it with the case as a bite relation, securing it to the arches with some blue mousse or similar low viscosity bite material applied when patient is in new accepted VDO.
- Written complete doctor's instructions to the lab.

*A 5% discount will be applied to the finished case if all the items are received correctly.*

If all information has been received correctly, the "in lab" time should be 2 weeks.