

MIS Seven	0.050" (1.27mm) Hex	1-72 UNF	Titanium	30 Ncm
Neoss				
Neoss	Neoss Proprietary	M2 X 0.40	Gold Alloy	32 Ncm
3.25mm	Neoss	M2 X 0.40	Gold Alloy	32 Ncm
NobelActive				
NobelActive 3.0	Unigrip	M1.4 X 0.30	Torqtite	15 Ncm
NobelActive NP	Unigrip	M1.6 X 0.35	Torqtite	35 Ncm
NobelActive RP/WP	Unigrip	M2X 0.40	Torqtite	35 Ncm
NobelReplace				
NobelReplace NP	Unigrip	M1.8 x 0.35	Titanium	35 Ncm
NobelReplace RP	Unigrip	M2 X 0.40	TorqTite	35 Ncm
OCO Biomedical	0.050" (1.27mm) Hex	1-72 UNF	Titanium	30 Ncm
Straumann Tissue Level				
Narrow Neck	SCS	M1.8 x 0.35	Titanium Alloy	35 Ncm
Regular Neck	SCS	M2 X 0.40	Titanium Alloy	35 Ncm
Wide Neck	SCS	M2 X 0.40	Titanium Alloy	35 Ncm
Straumann Bone-Level				
Narrow Cross-fit	SCS	M1.6 x 0.35	Titanium Alloy	35 Ncm
Regular Cross-fit	SCS	M1.6 x 0.35	Titanium Alloy	35 Ncm
Sybron Pro				
Sybron Pro (Internal Octa)	0.050" (1.27mm) Hex	M2 X 0.40	Titanium	35 Ncm
Sybron Pro (internal Hex)	0.050" (1.27mm) Hex	M2 X 0.40	Titanium	35 Ncm
Zimmer Spline	0.050" (1.27mm) Hex	1-72 UNF	Titanium	30 Ncm
Zimmer Screw-Vent	0.050" (1.27mm) Hex	1-72 UNF	Titanium	30 Ncm

Implant System	Prosthetic Screw Connection	Prosthetic Screw Thread	Prosthetic Screw Surface	Mfg. Recommended torque
Inclusive				30-35 Ncm
3I Certain Internal	0.048" (1.22mm) Hex	M1.6 X 0.35	Gold-Tite	20 Ncm
3I External	Square	M2 X 0.40	Gold-Tite	32-35 Ncm
Ankylos C/X	1mm Hex	M1.8 x 0.35	Titanium	15 Ncm
Astra EV				25 Ncm
Astra OsseoSpeed				
3.0 (Yellow)	0.050" (1.27mm) Hex	M1.4 X 0.30	Titanium	15 Ncm
3.5/4.0 (Aqua)	0.050" (1.27mm) Hex	M1.6 X 0.35	Titanium	20 Ncm
4.5/5.0 (Lilac)	0.050" (1.27mm) Hex	M2 X 0.40	Titanium	25 Ncm
BioHorizon				
BioHorizon Internal/External	0.050" (1.27mm) Hex		Titanium	30 Ncm
Branemark				
Branemark System NP	Unigrip	M1.6 X 0.35	Titanium	35 Ncm
Branemark System RP	Unigrip	M2 X 0.40	TorqTite	35 Ncm
Branemark System WP	Unigrip	M2.5 X 0.45	TorqTite	35 Ncm
Camlog ScrewLine				
3.3, 3.8, 4.3mm	0.050" (1.27mm) Hex	M1.6 x 0.35	Titanium	20 Ncm
5.0, 6.0mm	0.050" (1.27mm) Hex	M2 X 0.40	Titanium	20 Ncm
Endopore (Sybron)	0.050" (1.27mm) Hex		Titanium	30 Ncm
Hiossen				
Mini 3.5	1.2 mm Hex		EbonyGold, Titanium	20 Ncm
Standard 4.0, 4.5, 5.0, 6.0, 7.0	1.2 mm Hex		EbonyGold, Titanium	30 Ncm
Inclusive Tapered	0.050" (1.27mm) Hex	1-72 UNF		35 Ncm
Keystone Prima Connex	Square	M1.8 X 0.35	Gold	30 Ncm
Keystone Restore				
SD Restore	0.048" (1.22mm) Hex, Square	1-72 UNF	Titanium (Hex), Gold (Square)	30 Ncm
RD Restore	0.048" (1.22mm) Hex, Square	M2 x 0.40	Titanium (Hex), Gold (Square)	30 Ncm
WD Restore	0.048" (1.22mm) Hex, Square	M2.5 X 0.45	Titanium (Hex), Gold (Square)	30 Ncm
Keystone Renova	0.048" (1.22mm) Hex	M1.8 x 0.35	Titanium	30 Ncm
MIS Seven	0.050" (1.27mm) Hex	1-72 UNF	Titanium	30 Ncm
Neoss				
Neoss	Neoss Proprietary	M2 X 0.40	Gold Alloy	32 Ncm
3.25mm	Neoss	M2 X 0.40	Gold Alloy	32 Ncm
NobelActive				
NobelActive 3.0	Unigrip	M1.4 X 0.30	Torqtite	15 Ncm
NobelActive NP	Unigrip	M1.6 X 0.35	Torqtite	35 Ncm
NobelActive RP/WP	Unigrip	M2X 0.40	Torqtite	35 Ncm
NobelReplace				
NobelReplace NP	Unigrip	M1.8 x 0.35	Titanium	35 Ncm
NobelReplace RP	Unigrip	M2 X 0.40	TorqTite	35 Ncm